

## Liko Golvo 9000LB Lifter



- The Golvo LowBase (46 mm) allows for even greater clearance under low beds or stretchers versus the standard model
- Vertical lifting motion with flexible lift strap is ideal for horizontal lifting, ambulation and lifting from the floor
- Parallel base-width adjustment contributes to easier lifting of patients or residents to/from chairs, beds and the floor
- Armrests help to support patients during ambulation, as well as making it easier for caregivers to manoeuvre the lift
- User display highlights device usage, including lift count, service required indicator and overloads
- Optional quick-release hook system means caregivers can easily and rapidly change a sling bar or install a scale
- Ergonomic hand control for lift operation displays both battery usage and overload indicators
- An extensive assortment of accessories is available including slings, slingbars and scales

### Warranty Information

- Frames: 10 Years
- Electrical Parts: 2 Years

Contract approved

### Support @ Home

**Function:** Lifting persons

**Category:** Prescribed

Click [here](#) for more information

### NDIS

**Item Code:**

05\_123603161\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material
LSS393105	Liko Golvo 9000LB Lifter	1185 mm	745 - 1028 mm	1455 - 2090 mm	200 kg	Aluminium

