

## BodyPoint Calf Strap

- Aeromesh TM fabric allows for air circulation and provides comfortable support strap for legs
- Designed to attach to legrest hangers of wheelchairs for lower leg support

 Contract approved

### Support @ Home

**Function:** Wheelchair accessories

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- Limited Lifetime



Code	Display Name	Total Length
WCA860464	Calf Strap - Small	800 mm
WCA860465	Calf Strap - Medium	900 mm
WCA860466	Calf Strap - Large	1000 mm