

Aspire Toilet Seat Raiser



- Assists toileting and transfers by raising the toilet seat height
- Moulded plastic seat sits securely over the existing toilet
- Fits all standard toilet seats
- Adjustable mounting brackets to allow for correct fitting to different size toilets
- Available in three different height versions to suit individual requirements
- SWL: 200 kg

Warranty Information

- 1 Year



✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_091218063_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information

Compliance

AS/NZS 3973:2009

| Code | Display Name | RRP | Box Quantity |
|-----------|-----------------------------------|------|--------------|
| BTT149800 | Aspire Toilet Seat Raiser - 50mm | 69.5 | 1 |
| BTT150200 | Aspire Toilet Seat Raiser - 100mm | 75 | 1 |
| BTT150300 | Aspire Toilet Seat Raiser - 150mm | 82 | 1 |