

Aspire Toilet Seat Raiser



- Assists toileting and transfers by raising the toilet seat height
- Moulded plastic seat sits securely over the existing toilet
- Fits all standard toilet seats
- Adjustable mounting brackets to allow for correct fitting to different size toilets
- Available in three different height versions to suit individual requirements
- SWL: 200 kg

Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

ndis

Item Code:

05_091218063_0103_1_2

Toilet Aids – Composite:

05_414_0103_1_2

Click [here](#) for more information



Compliance

AS/NZS 3973:2009

Code	Display Name	Material	Unit Weight
BTT150300	Aspire Toilet Seat Raiser - 150mm	Plastic	1.4 kg