

Invacare Birdie Compact Patient Lifter



- Birdie “Compact” has a 55 cm width and can be folded and unfolded without the need for tools or special skills required. This makes it ideal for the nursing or home care environments where space and time are at a premium
- The Birdie Compact is the only lifter that is TÜV certified that is this narrow – unique to the market
- All attachments on the Birdie Compact are quick release, making it very easy to change when required
- Intelligent control system includes a service reminder, emergency operation, audio battery alarm, and battery indicator
- Manual leg width adjustment
- Safe patient handling with up to 150 kg
- Supplied with yoke attachment
- Electric emergency lowering up and down
- Item Weight: 36 kg

Optional Accessories:

- Pivot Frame (LSS364112)
- Full range of slings

Warranty Information

- 3 Years Electronics
- 2 Years Battery
- 1 Year Handset

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range	Leg Gait Range
LSS364100	Invacare Birdie Compact Patient Lifter	1000 mm	550 mm	1600 mm	150 kg	Steel	36 kg	530 - 1600 mm	550 - 890 mm