

DermaSaver Stay-Put Heel Protector



- The DermaSaver Stay-Put Heel Protector protects fragile skin from pressure, friction and minor traumas that result in tears, bruising and minor abrasion
- The heel protector is soft and loose with ample air-flow reaching the skin
- It is constructed of two layers of soft, lightweight microspring textile, which allows air-flow and moisture transfer
- Anti-microbial - will help prevent growth of bacteria, latex free and fully washable
- Measure ankle on a diagonal from heel to top of foot. For borderline measurements, choose the larger size for extra comfort

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Ankle Circumference
PTA533313	DermaSaver Stay-Put Heel Protector - XS	< 280 mm
PTA533313S	DermaSaver Stay-Put Heel Protector - S	280 - 300 mm
PTA533313M	DermaSaver Stay-Put Heel Protector - M	300 - 370 mm
PTA533313L	DermaSaver Stay-Put Heel Protector - L	370 - 390 mm
PTA533313XL	DermaSaver Stay-Put Heel Protector - XL	> 390 mm