

MOTOmed Gracile 12



RECK
MOTOmed®

Effective movement therapy for children with special needs, it is the only motorised movement therapy device that fulfills the special needs of children.



- The Gracile12 comes with both upper limb and lower limb trainers installed
- During active training, children can play MOTOMax for more fun during training.
- Leg trainer (optional: arm/upper body trainer).
- Paediatric safety foot shells with paediatric leg guides.
- Small distance of 12 cm between the inner rims of the foot shells.
- Two-stage pedal radius adjustment (3.5 cm / 7.0 cm).
- Height adjustable pedal axle.
- Handlebar adjustable in height and inclination.
- Mobile with big castors (13cm diameter).
- All-metal construction, premium and stable.
- Additional stability by telescopic front stand.
- MOTOMed SmoothDriveSystem with a gentle startup and rundown.
- MovementProtector & SpasmControl.
- Speed regulation from 0 to 60 rpm (passive).
- Regulation of motor power from 1 to 10.
- Preset the therapy duration from 0 to 120 minutes.
- User-friendly operating panel with a colour coded system and big, easy-feel buttons.
- Language selection between 27 different languages.

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight
FPM101112	MOTomed Gracile 12 - Custom (as per script)	60 - 88 cm	55 cm	99 - 114 cm	95 kg	33 kg
FPM101110	MOTomed Gracile w/ Arms - Custom (as per script)	600 - 880 mm	550 mm	990 - 1140 mm	90 Kg	33 Kg