

## Conni Reusable Bed Pad with Tuck-ins



*Conni.*

- The Conni absorbent and waterproof bed pads with tuck-ins assist in the management of incontinence, bedwetting and toilet training
- Designed to hold 2.5 litres over an eight hour period, which means they are the perfect addition to any single or king-single beds that require full void protection
- Cost-effective and environmentally sustainable alternative to disposable.
- Certified as being free from harmful materials and are safe even against newborn skin
- Soft and breathable
- Washable and reusable

### Warranty Information

- 1 Year



**ndis**

**Item Code:**

03\_092489060\_0103\_1\_1

Click [here](#) for more information

**Compliance**

No applicable compliance standard

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Colour	Capacity
CON228510	Conni Bed Pad with Tuck-ins - Mauve	81	11	1000 mm	1000 mm	Mauve	2.5 litres
CON228515	Conni Bed Pad with Tuck-ins - Teal Blue	81	11	1000 mm	1000 mm	Teal Blue	2.5 litres
CON228516	Conni Bed Pad with Tuck-ins - White	81	-	1000 mm	1000 mm	White	2.5 litres
CON228527	Conni Bed Pad with Tuck-ins - Animal print	81	-	1000 mm	1000 mm	Animal print	2.5 litres