

TheraBand Tubing Retail Pack



- Ideal for upper body & lower body Exercise to strengthen muscles, increase range of motion, and improve flexibility.
- Convenient and portable product that works in the clinic, the gym, at home or on-the-go.
- Color-coded Progress
- x3 in a pack:
 - Light pack includes: Low, Medium, Heavy
 - Heavy pack includes: Heavy, Extra Heavy, Special Heavy



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Colour	Resistance
FPP311100	TheraBand Tubing Retail Pack - Light	1.83 m	140 mm	-	Heavy
FPP311200	TheraBand Tubing Retail Pack - Heavy	1.83 m	140 mm	Multi	Heavy