

## Elephant Feet Bed/Chair Raiser



- Quality plastic raiser used to increase the height of a bed or chair
- Ideal for persons who have difficulty standing up from a bed or chair
- Dished recess for added stability
- Supplied with non-slip rubber insert
- Sold individually
- Recess diameter 90 mm
- Maximum load 380 kg (per set)

✔ Contract approved

### Support @ Home

**Function:** Height adjustment of furniture

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Total Height
BEA011701	Elephant Feet Bed/Chair Raiser - 140mm	140 mm
BEA011705	Elephant Feet Bed/Chair Raiser - 140mm - Set of 4	140 mm