

Elephant Feet Bed/Chair Raiser



- Quality plastic raiser used to increase the height of a bed or chair
- Ideal for persons who have difficulty standing up from a bed or chair
- · Dished recess for added stability
- Supplied with non-slip rubber insert
- Sold individually
- Recess diameter 90 mm
- Maximum load 380 kg (per set)

Display Name	Total Height
Elephant Feet Bed/Chair Raiser - 90mm	90 mm
Elephant Feet Bed/Chair Raiser - 140mm	140 mm
	lephant Feet Bed/Chair Raiser - 90mm

BEA011705 Elephant Feet Bed/Chair Raiser - 140mm - Set of 4 140 mm

Ocontracted DVA produ	ct

Ocontract approved