

Anaesthetic Stool



- Comfortable, height adjustable stool for use in the workplace
- 50 mm twin wheels
- Gas assisted height adjustment
- Padded vinyl seat
- Polished aluminium five star base
- Adjustable height: 49 cm - 63 cm
- Cushion diameter: 35 cm
- Foam thickness: 6 cm
- Premium version has additional foot ring



Compliance

No applicable compliance standard

Warranty Information

- 3 Years

Code	Display Name
WAT803100	Anaesthetic Stool Premium with Foot Ring - Navy
WAT803200	Anaesthetic Stool Standard - Navy