

Compreflex No Foot & Basic Liners



- The Compreflex Standard Calf is an inelastic compression device. Paired with the non-compressive Basic Liner, the Compreflex Standard Calf offers compression to the calf.
- Accurate compression measurements of 20-30, 30-40, 40-50 powered by the patent pending.
- AccutabEasy-to-use features including the built in stretch panel that holds the product in place before straps are closed.
- Product contains a pair of non-compressive underliners for the protection of the skin.
- Product material is short stretch, allowing for a firm compressive hold while calf contracts and expands.
- Durable material that allows for long-term wear.
- The Compreflex Standard Calf is ideal for patients with Venous and Lymphatic disorders.
- Patients with Chronic Venous Insufficiency resulting in venous leg ulcers or those with chronic edema will find the product easy to use and supportive in their care.
- Individuals with lymphedema, actively working to manage their condition will find the conforming material and supporting application features of the Compreflex ideal for the management of the disorder.

✓ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Colour	Calf/Arm Circumference	Ankle Circumference
FBS310200	Compreflex No Foot & Basic Liners - Large - Tall Black	340 - 420 mm	Black	400 - 500 mm	260 - 360 mm
FBS310210	Compreflex No Foot & Basic Liners - Medium - Reg Black	300 - 339 mm	Black	340 - 440 mm	210 - 300 mm

Code	Display Name	Total Length	Colour	Calf/Arm Circumference	Ankle Circumference
FBS310220	Compreflex No Foot & Basic Liners - Medium - Tall Black	340 - 420 mm	Black	340 - 440 mm	210 - 300 mm
FBS310230	Compreflex No Foot & Basic Liners - XL- Tall Black	340 - 420 mm	Black	310 - 410 mm	310 - 410 mm