

Aspire ComfiMotion Fixed Companion Bed




- Matching single bed without any adjustability designed to sit flush when positioned together with the ComfiMotion CARE Bed and ComfiMotion Activ Care Bed
- It facilitates a side-by-side sleeping arrangement at a fixed height of 370 mm, without compromising on the users care and assistance requirements

Warranty Information

- 10 years Frame
- 1 year Fabric

Ocontracted		
Contracted	DVA	product



Compliance

IEC 60601-2-52:2009



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEB046800 A	Aspire ComfiMotion Fixed Companion Bed - Portobello Grey - Long Single (Headboard Se	p) 2030 mm	920 mm	370 mm	180 kg