

Aspire ComfiMotion Activ Care Bed Package - Chesterfield



The Aspire ComfiMotion Activ Care Bed Package is designed to deliver advanced support, flexibility and style in one complete solution. Supplied with the Chesterfield Free Standing Headboard as standard, this package combines the clinical functionality of the Activ Care Bed with a modern, homely look - supporting both quality care and a welcoming bedroom environment.

The Activ Care bed combines the best of in-home care functionality, positioning and low height safety with luxurious and contemporary style for enhanced comfort, care and wellbeing.

- **Ultra Low Height:** Minimum platform height of 195 mm reduces falls risk, improving user safety
- **Customisable Safe Entry/Exit Height:** Programmable single button press ensures perfect bed height to maximise independence
- **Under Bed Clearance:** provides clear access for patient lifters and overbed tables
- **Battery Backup:** Provides confidence and maintains independence during power outages
- **Night Assist:** Luminous underbed halo for safe entry and exit
- **Illuminated Handset:** Large button interface is intuitive and highly visible at all times
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief
- **Hi-Lo Functionality:** Assists with manual handling and helps to mitigate fall risks
- **Head Tilt:** Relaxes neck muscles, providing optimal support whilst reading or watching TV
- **Back Raise:** Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea

 Contract approved

- **Auto-Regression in the Backrest:** Stops the user sliding down the bed during movement while improving comfort and maximising pressure redistribution
- **Clinical Positioning:** One touch raises back and legs to alleviate cardiac exertion, improving resting comfort
- **Leg Raise:** Promotes circulation, relieving tired, swollen/ fatigued or painful limbs
- **Memory Function:** Allows the user to set and return to their preferred positions at any time with one click of a button
- **Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back
- **Seated Position:** Assists with bed entry and exit and upright dining
- **Trendelenburg (and Reverse):** Can assist care givers when repositioning a user
- **Lockable Castors:** Eight heavy duty dual castors
- Accessories sold separately

*Carers should seek advice from a trained Allied Health professional for correct manual handling techniques.

Warranty Information

- 10 Years Frame
- 2 Years Electronics
- 1 Year Handset, Cables & Fabric
- 1 Year Headboard

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Colour
BEB046770V2	Aspire ComfiMotion Activ Care Bed Package - Chesterfield - Queen (inc. Freestanding Headboard)	2030 mm	1520 mm	195 - 700 mm	360 kg	156 kg	Portobello Grey
BEB046740V2	Aspire ComfiMotion Activ Care Bed Package - Chesterfield - Long Single (inc. Freestanding Headboard)	2030 mm	920 mm	195 - 700 mm	250 kg	114 kg	Portobello Grey
BEB046750V2	Aspire ComfiMotion Activ Care Bed Package - Chesterfield - King Single (inc. Freestanding Headboard)	2030 mm	1070 mm	195 - 700 mm	250 kg	126 kg	Portobello Grey
BEB046760V2	Aspire ComfiMotion Activ Care Bed Package - Chesterfield - Long Double (inc. Freestanding Headboard)	2030 mm	1370 mm	195 - 700 mm	360 kg	149 kg	Portobello Grey