

Penny Wingback Chair



- Modern and stylish design with a contemporary take on a classic wingback chair
- Versatile design that will complement any interior setting
- Pair this armchair with a side table for an inviting place to sit and read
- Narrow solid timber legs with a walnut finish provide a sturdy base
- Quality construction ensures this armchair will withstand daily use
- Fabric colour/material to be advised at time of order



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
FUL052180	Penny Wingback Chair	770mm	660mm	1105mm