

Aspire Deluxe General Purpose Sling with Head Support - Mesh - Black Clips



- Clip fitting option to fit sling to a pivot frame
- Australian standards tested and registered with the TGA
- Versatile sling offering body contoured fabric shaping for maximum patient comfort and support
- Durable keyhole clip attachment points for safety and security during transfers
- This sling has a comfortable padded head support with rigidisers for posture management
- Quality - manufactured in the UK using European materials, computer guided stitching and CNC fabric cutting
- User Comfort & Dignity - features additional padding in weight bearing areas and option to position strap around the outside of legs for more discreet positioning (applicable models)
- Serial Numbered - slings are embroidered with a unique serial number to maintain traceability and compliance at all times
- Service & Maintenance - service record labels ensure you can quickly and easily identify service history
- Easy Sling Identification - all slings are colour coded by size and have reinforced edging for longevity
- Easy Clean - machine washable up to 85°C temperature
- For FIM score 1 or 2

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 1 Year

Code	Display Name	Safe Working Load	Material
LSS390916B	Deluxe General Purpose Sling with Head Support - Mesh - Black Clips - Small	200 kg	Mesh
LSS390917B	Deluxe General Purpose Sling with Head Support - Mesh - Black Clips - Medium	200 kg	Mesh
LSS390918B	Deluxe General Purpose Sling with Head Support - Mesh - Black Clips - Large	300 kg	Mesh
LSS390919B	Deluxe General Purpose Sling with Head Support - Mesh - Black Clips - Extra Large	300 kg	Mesh