

Aspire Autoclavable Kidney Dishes




Our range of Kidney Dishes are fully autoclavable and made from premium polypropylene. They will last for many years in any healthcare environment subject to the rigours of cleaning, sanitation and infection control requirements.

- Heavy duty and reusable
- Designed for use in healthcare environments such as surgeries, acute care, general ward, medical centres, and aged care facilities
- Autoclavable polypropylene that can be cleaned using autoclaves and sanitisers
- Pastel green colour for easy identification

Warranty Information

- 1 Year



 Contract approved

Compliance

No applicable compliance standard

Code	Display Name	Total Length	Capacity
WAT770400	Kidney Dish Clear Pastic Disposable 210mm (700ml) -		700 ml
WAT770500	Kidney Dish - Small	160 mm	200 ml
WAT763700	Kidney Dish - Medium	220 mm	400 ml
WAT770600	Kidney Dish - Large	255 mm	600 ml
WAT770700	Kidney Dish - Extra Large	300 mm	800 ml