

MEDI K.620 PROTECT GENU ELASTIC KNEE SUPPORT WITH PATELLA RING



- Silicone patellar ring for massaging and reducing swelling.
- Rods incorporated into both sides of the support for increased stabilization.
- Increased comfort by using softer material in the flexion area.
- Breathable, elastic and light material, and anatomically shaped knit.

Application:

- Massage effect due to silicone insert
- Improves the blood circulation and aids absorption of haematomas and oedemas
- Improves proprioception and muscular support (stabilization of the knee)

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS307340	MEDI K.620 PROTECT GENU ELATIC KNEE SUPPORT WITH PATELLA RING - SIZE II
FBS307350	MEDI K.620 PROTECT GENU ELATIC KNEE SUPPORT WITH PATELLA RING - SIZE III
FBS307360	MEDI K.620 PROTECT GENU ELATIC KNEE SUPPORT WITH PATELLA RING - SIZE IV
FBS307370	MEDI K.620 PROTECT GENU ELATIC KNEE SUPPORT WITH PATELLA RING - SIZE V
FBS307380	MEDI K.620 PROTECT GENU ELATIC KNEE SUPPORT WITH PATELLA RING - SIZE VI
FBS307390	MEDI K.620 PROTECT GENU ELATIC KNEE SUPPORT WITH PATELLA RING - SIZE VII