

Stocking Aid - JOBST Stocking Donner



- A plastic coated metal stocking applicator to assist in putting on compression or surgical stockings
- The frame has foam covered handles
- The frame will hold the stocking/sock sufficiently open to allow the foot to be inserted
- The stocking still needs to be fed over the frame using the hands but may eliminate prolonged bending to reach the feet
- Available in one size and colour

Warranty Information

• 1 Year



Ocontract approved

ndis

Item Code:

05_090900046_0103_1_2 Click **here** for more information

Compliance

No applicable compliance standard



Code Display Name

PTS641090 Stocking Aid - JOBST Stocking Donner