

Stocking Aid - JOBST Stocking Donner



- A plastic coated metal stocking applicator to assist in putting on compression or surgical stockings.
- The frame has foam covered handles.
- The frame will hold the stocking/sock sufficiently open to allow the foot to be inserted.
- The stocking still needs to be fed over the frame using the hands but may eliminate prolonged bending to reach the feet.
- Available in one size and colour.

Warranty Information

- 1 Year



Code	Display Name
PTS641090	Stocking Aid - JOBST Stocking Donner

Contracted DVA product

Contract approved

ndis

Item Code:

05_090900046_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard