

## **TabTimer Rainbow CONNECTION Pill Box Timer**

## **TabTimer Rainbow CONNECTION Pill Box Timer**





- Up-to 4 alarms per day
- 1-15 compartments for storing different medications
- Easy to set and use
- When due, alarms sounds for 1 minute
- After 5 minutes the missed pill reminder sounds for 30 seconds a total of 3 times
- 12 / 24 hour clock
- LCD display
- Volume level approx 45-50db
- 2 x AAA size battery included
- Food Grade Tritan materials BPA-free
- Weight 83g
- Timer Only Dimensions: 70mm x 70mm x 25mm
- Timer with 7 Days connected Dimensions: 70mm x 300mm x 25mm
- Each compartment is approx. 25mm x 25mm x 22 mm (large enough for approx. 6-8 large vitamin tablets, 16 panadol, 10-25+ smaller pills)
- 12 Month Limited Warranty by product source (excludes battery)

Code

**Display Name** 

DLG275050 TabTimer Rainbow CONNECTION Pill Box Timer