

Aspire ComfiMotion Bed Collection - Accessibility Castor (each)



- Larger castors raise the bed frame to a minimum clearance of 155 mm facilitating patient lifter access
- Designed to be used in conjunction with ComfiMotion bed collection
- x4 required for a complete bed

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Not on S@H

Category: Not on S@H

Click [here](#) for more information



Code	Display Name
MWT451080	Aspire ComfiMotion Bed Collection - Accessibility Castor (each)
GSP640100	Low Profile Fixed Leg - Aspire ComfiMotion Bed Collection