

Aspire ComfiMotion Bed Collection - Accessibility Castor (each)



- Larger castors raise the bed frame to a minimum clearance of 155 mm facilitating patient lifter access
- Designed to be used in conjunction with ComfiMotion bed collection
- x4 required for a complete bes

- ✓ Contracted DVA product
- ✓ Contract approved



Code	Display Name	RRP	Box Quantity
MWT451080	Aspire ComfiMotion Bed Collection - Accessibility Castor (each)	50.5	0
GSP640100	Low Profile Fixed Leg - Aspire ComfiMotion Bed Collection	43.8	-