

Sigvaris Cotton Compression Stocking



- Unisex compression stockings for people with sensitive skin
- The mix of fibres also ensures that the stocking slides on and off easily
- Excellent thermoregulation
- Opaque with a subtle ribbed texture
- Tangible high quality as the result of using precious natural fibres
- Virtually lint-free cotton fabric with good stretching properties
- Compression: CCL 1 (18 - 21 mmHg), CCL 2 (23 - 32 mmHg), CCL 3 (34 - 46 mmHg), CCL 4 (=49 mmHg)
- Please specify at the time of ordering:
 - Size: XS, S, M, L (Please refer to images for size guide)
 - Toe: Open or Closed

✓ Contracted DVA product

✓ Contract approved

Compliance

No applicable compliance standard

Code	Display Name	Colour
PTS648110	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Short Nature Open Toe	Natural
PTS648150	SIGVARIS Cotton - Knee - 23-32mmHg XXS Norm Nature Open Toe	Natural
PTS648145	SIGVARIS Cotton - Knee - 23-32mmHg XXS Norm Nature Closed Toe	Natural
PTS648140	SIGVARIS Cotton - Knee - 23-32mmHg XXS Norm Black Open Toe	Black
PTS648135	SIGVARIS Cotton - Knee - 23-32mmHg XXS Norm Black Closed Toe	Black
PTS648130	SIGVARIS Cotton - Knee - 23-32mmHg XXS Long Nature Open Toe	Natural
PTS648125	SIGVARIS Cotton - Knee - 23-32mmHg XXS Long Nature Closed Toe	Natural

Code	Display Name	Colour
PTS648065	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Long Nature Closed Toe	Natural
PTS648115	SIGVARIS Cotton - Knee - 23-32mmHg XXS Long Black Closed Toe	Black
PTS648165	SIGVARIS Cotton - Knee - 23-32mmHg XXS Short Nature Closed Toe	Natural
PTS648105	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Short Nature Closed Toe	Natural
PTS648100	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Short Black Open Toe	Black
PTS648095	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Short Black Closed Toe	Black
PTS648090	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Norm Nature Open Toe	Natural
PTS648085	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Norm Nature Closed Toe	Natural
PTS648080	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Norm Black Open Toe	Black
PTS648075	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Norm Black Closed Toe	Black
PTS648680	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg L Short Nature Open Toe	Natural
PTS648120	SIGVARIS Cotton - Knee - 23-32mmHg XXS Long Black Open Toe	Black
PTS648205	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Norm Nature Closed Toe	Natural
PTS648255	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg M PLUS Long Nature Closed Toe	Natural
PTS648250	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg L PLUS Norm Nature Open Toe	Natural
PTS648245	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg L PLUS Norm Nature Closed Toe	Natural
PTS648240	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg L PLUS Long Nature Open Toe	Natural

Code	Display Name	Colour
PTS648235	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg L PLUS Long Nature Closed Toe	Natural
PTS648230	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Short Nature Open Toe	Natural
PTS648225	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Short Nature Closed Toe	Natural
PTS648220	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Short Black Open Toe	Black
PTS648155	SIGVARIS Cotton - Knee - 23-32mmHg XXS Short Black Closed Toe	Black
PTS648210	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Norm Nature Open Toe	Natural
PTS648160	SIGVARIS Cotton - Knee - 23-32mmHg XXS Short Black Open Toe	Black
PTS648200	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Norm Black Open Toe	Black
PTS648195	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Norm Black Closed Toe	Black
PTS648190	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Long Nature Open Toe	Natural
PTS648185	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Long Nature Closed Toe	Natural
PTS648180	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Long Black Open Toe	Black
PTS648175	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Long Black Closed Toe	Black
PTS648170	SIGVARIS Cotton - Knee - 23-32mmHg XXS Short Nature Open Toe	Natural
PTS648060	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Long Black Open Toe	Black
PTS648215	SIGVARIS Cotton - Knee - 23-32mmHg XXS PLUS Short Black Closed Toe	Black
PTS647905	SIGVARIS Cotton - Knee - 23-32mmHg XL Norm Nature Closed Toe	Natural
PTS647950	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Norm Nature Open Toe	Natural

Code	Display Name	Colour
PTS647945	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Norm Nature Closed Toe	Natural
PTS647940	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Norm Black Open Toe	Black
PTS647935	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Norm Black Closed Toe	Black
PTS647930	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Long Nature Open Toe	Natural
PTS647925	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Long Nature Closed Toe	Natural
PTS647920	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Long Black Open Toe	Black
PTS648070	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Long Nature Open Toe	Natural
PTS647910	SIGVARIS Cotton - Knee - 23-32mmHg XL Norm Nature Open Toe	Natural
PTS647965	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Short Nature Closed Toe	Natural
PTS647900	SIGVARIS Cotton - Knee - 23-32mmHg XL Norm Black Open Toe	Black
PTS647895	SIGVARIS Cotton - Knee - 23-32mmHg XL Norm Black Closed Toe	Black
PTS647890	SIGVARIS Cotton - Knee - 23-32mmHg XL Long Nature Open Toe	Natural
PTS647885	SIGVARIS Cotton - Knee - 23-32mmHg XL Long Nature Closed Toe	Natural
PTS647880	SIGVARIS Cotton - Knee - 23-32mmHg XL Long Black Open Toe	Black
PTS647875	SIGVARIS Cotton - Knee - 23-32mmHg XL Long Black Closed Toe	Black
PTS647870	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Short Nature Open Toe	Natural
PTS647865	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Short Nature Closed Toe	Natural
PTS647915	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Long Black Closed Toe	Black

Code	Display Name	Colour
PTS648005	SIGVARIS Cotton - Knee - 23-32mmHg XS Long Nature Closed Toe	Natural
PTS648055	SIGVARIS Cotton - Knee - 23-32mmHg XS PLUS Long Black Closed Toe	Black
PTS648050	SIGVARIS Cotton - Knee - 23-32mmHg XS Short Nature Open Toe	Natural
PTS648045	SIGVARIS Cotton - Knee - 23-32mmHg XS Short Nature Closed Toe	Natural
PTS648040	SIGVARIS Cotton - Knee - 23-32mmHg XS Short Black Open Toe	Black
PTS648035	SIGVARIS Cotton - Knee - 23-32mmHg XS Short Black Closed Toe	Black
PTS648030	SIGVARIS Cotton - Knee - 23-32mmHg XS Norm Nature Open Toe	Natural
PTS648025	SIGVARIS Cotton - Knee - 23-32mmHg XS Norm Nature Closed Toe	Natural
PTS648020	SIGVARIS Cotton - Knee - 23-32mmHg XS Norm Black Open Toe	Black
PTS647955	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Short Black Closed Toe	Black
PTS648010	SIGVARIS Cotton - Knee - 23-32mmHg XS Long Nature Open Toe	Natural
PTS647960	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Short Black Open Toe	Black
PTS648000	SIGVARIS Cotton - Knee - 23-32mmHg XS Long Black Open Toe	Black
PTS647995	SIGVARIS Cotton - Knee - 23-32mmHg XS Long Black Closed Toe	Black
PTS647990	SIGVARIS Cotton - Knee - 23-32mmHg XL Short Nature Open Toe	Natural
PTS647985	SIGVARIS Cotton - Knee - 23-32mmHg XL Short Nature Closed Toe	Natural
PTS647980	SIGVARIS Cotton - Knee - 23-32mmHg XL Short Black Open Toe	Black
PTS647975	SIGVARIS Cotton - Knee - 23-32mmHg XL Short Black Closed Toe	Black

Code	Display Name	Colour
PTS647970	SIGVARIS Cotton - Knee - 23-32mmHg XL PLUS Short Nature Open Toe	Natural
PTS648270	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg M PLUS Norm Nature Open Toe	Natural
PTS648015	SIGVARIS Cotton - Knee - 23-32mmHg XS Norm Black Closed Toe	Black
PTS648520	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg M Norm Nature Open Toe	Natural
PTS648565	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XL Long Nature Closed Toe	Natural
PTS648560	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg S Short Nature Open Toe	Natural
PTS648555	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg S Short Nature Closed Toe	Natural
PTS648550	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg S Norm Nature Open Toe	Natural
PTS648545	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg S Norm Nature Closed Toe	Natural
PTS648540	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg S Long Nature Open Toe	Natural
PTS648535	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg S Long Nature Closed Toe	Natural
PTS648260	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg M PLUS Long Nature Open Toe	Natural
PTS648525	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg M Short Nature Closed Toe	Natural
PTS648580	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XL Norm Nature Open Toe	Natural
PTS648515	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg M Norm Nature Closed Toe	Natural
PTS648510	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg M Long Nature Open Toe	Natural
PTS648505	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg M Long Nature Closed Toe	Natural
PTS648500	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg L Short Nature Open Toe	Natural

Code	Display Name	Colour
PTS648495	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg L Short Nature Closed Toe	Natural
PTS648490	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg L Norm Nature Open Toe	Natural
PTS648485	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg L Norm Nature Closed Toe	Natural
PTS648480	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg L Long Nature Open Toe	Natural
PTS648530	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg M Short Nature Open Toe	Natural
PTS648620	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XS Short Nature Open Toe	Natural
PTS648670	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg L Norm Nature Open Toe	Natural
PTS648665	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg L Norm Nature Closed Toe	Natural
PTS648660	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg L Long Nature Open Toe	Natural
PTS648655	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg L Long Nature Closed Toe	Natural
PTS648650	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XXS Short Nature Open Toe	Natural
PTS648645	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XXS Short Nature Closed Toe	Natural
PTS648640	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XXS Norm Nature Open Toe	Natural
PTS648635	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XXS Norm Nature Closed Toe	Natural
PTS648570	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XL Long Nature Open Toe	Natural
PTS648625	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XXS Long Nature Closed Toe	Natural
PTS648575	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XL Norm Nature Closed Toe	Natural
PTS648615	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XS Short Nature Closed Toe	Natural

Code	Display Name	Colour
PTS648610	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XS Norm Nature Open Toe	Natural
PTS648605	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XS Norm Nature Closed Toe	Natural
PTS648600	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XS Long Nature Open Toe	Natural
PTS648595	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XS Long Nature Closed Toe	Natural
PTS648590	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XL Short Nature Open Toe	Natural
PTS648585	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XL Short Nature Closed Toe	Natural
PTS648465	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XXS Norm Nature Closed Toe	Natural
PTS648630	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg XXS Long Nature Open Toe	Natural
PTS648310	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XL PLUS Norm Nature Open Toe	Natural
PTS648475	SIGVARIS Cotton - Thigh w/Grip Top - 18-21mmHg L Long Nature Closed Toe	Natural
PTS648355	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg L Long Nature Closed Toe	Natural
PTS648350	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XXS PLUS Norm Nature Open Toe	Natural
PTS648345	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XXS PLUS Norm Nature Closed Toe	Natural
PTS648340	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XXS PLUS Long Nature Open Toe	Natural
PTS648335	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XXS PLUS Long Nature Closed Toe	Natural
PTS648330	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XS PLUS Norm Nature Open Toe	Natural
PTS648325	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XS PLUS Norm Nature Closed Toe	Natural
PTS648365	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg L Norm Nature Closed Toe	Natural

Code	Display Name	Colour
PTS648315	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XS PLUS Long Nature Closed Toe	Natural
PTS648370	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg L Norm Nature Open Toe	Natural
PTS648305	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XL PLUS Norm Nature Closed Toe	Natural
PTS648300	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XL PLUS Long Nature Open Toe	Natural
PTS648295	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XL PLUS Long Nature Closed Toe	Natural
PTS648290	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg S PLUS Norm Nature Open Toe	Natural
PTS648285	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg S PLUS Norm Nature Closed Toe	Natural
PTS648280	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg S PLUS Long Nature Open Toe	Natural
PTS648275	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg S PLUS Long Nature Closed Toe	Natural
PTS647850	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Norm Nature Open Toe	Natural
PTS648320	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XS PLUS Long Nature Open Toe	Natural
PTS648415	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XL Long Nature Closed Toe	Natural
PTS648265	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg M PLUS Norm Nature Closed Toe	Natural
PTS648460	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XXS Long Nature Open Toe	Natural
PTS648455	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XXS Long Nature Closed Toe	Natural
PTS648450	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XS Norm Nature Open Toe	Natural
PTS648445	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XS Norm Nature Closed Toe	Natural
PTS648440	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XS Long Nature Open Toe	Natural

Code	Display Name	Colour
PTS648435	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XS Long Nature Closed Toe	Natural
PTS648430	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XL Norm Nature Open Toe	Natural
PTS648360	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg L Long Nature Open Toe	Natural
PTS648420	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XL Long Nature Open Toe	Natural
PTS648470	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XXS Norm Nature Open Toe	Natural
PTS648410	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg S Norm Nature Open Toe	Natural
PTS648405	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg S Norm Nature Closed Toe	Natural
PTS648400	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg S Long Nature Open Toe	Natural
PTS648395	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg S Long Nature Closed Toe	Natural
PTS648390	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg M Norm Nature Open Toe	Natural
PTS648385	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg M Norm Nature Closed Toe	Natural
PTS648380	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg M Long Nature Open Toe	Natural
PTS648375	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg M Long Nature Closed Toe	Natural
PTS648425	SIGVARIS Cotton - Knee w/Grip Top - 23-32mmHg XL Norm Nature Closed Toe	Natural
PTS647285	SIGVARIS Cotton - Knee - 18-21mmHg L Long Nature Closed Toe	Natural
PTS647330	SIGVARIS Cotton - Knee - 18-21mmHg M Long Nature Open Toe	Natural
PTS647325	SIGVARIS Cotton - Knee - 18-21mmHg M Long Nature Closed Toe	Natural
PTS647320	SIGVARIS Cotton - Knee - 18-21mmHg M Long Black Open Toe	Black

Code	Display Name	Colour
PTS647315	SIGVARIS Cotton - Knee - 18-21mmHg M Long Black Closed Toe	Black
PTS647310	SIGVARIS Cotton - Knee - 18-21mmHg L Norm Nature Open Toe	Natural
PTS647305	SIGVARIS Cotton - Knee - 18-21mmHg L Norm Nature Closed Toe	Natural
PTS647300	SIGVARIS Cotton - Knee - 18-21mmHg L Norm Black Open Toe	Black
PTS647440	SIGVARIS Cotton - Knee - 18-21mmHg XS Long Black Open Toe	Black
PTS647290	SIGVARIS Cotton - Knee - 18-21mmHg L Long Nature Open Toe	Natural
PTS647345	SIGVARIS Cotton - Knee - 18-21mmHg M Norm Nature Closed Toe	Natural
PTS647280	SIGVARIS Cotton - Knee - 18-21mmHg L Long Black Open Toe	Black
PTS647275	SIGVARIS Cotton - Knee - 18-21mmHg L Long Black Closed Toe	Black
PTS647270	SIGVARIS Cotton - Knee - 18-21mmHg XXS PLUS Norm Nature Open Toe	Natural
PTS647265	SIGVARIS Cotton - Knee - 18-21mmHg XXS PLUS Norm Nature Closed Toe	Natural
PTS647260	SIGVARIS Cotton - Knee - 18-21mmHg XXS PLUS Norm Black Open Toe	Black
PTS647255	SIGVARIS Cotton - Knee - 18-21mmHg XXS PLUS Norm Black Closed Toe	Black
PTS647250	SIGVARIS Cotton - Knee - 18-21mmHg XXS PLUS Long Nature Open Toe	Natural
PTS647245	SIGVARIS Cotton - Knee - 18-21mmHg XXS PLUS Long Nature Closed Toe	Natural
PTS647295	SIGVARIS Cotton - Knee - 18-21mmHg L Norm Black Closed Toe	Black
PTS647385	SIGVARIS Cotton - Knee - 18-21mmHg S Norm Nature Closed Toe	Natural
PTS647860	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Short Black Open Toe	Black

Code	Display Name	Colour
PTS647430	SIGVARIS Cotton - Knee - 18-21mmHg XL Norm Nature Open Toe	Natural
PTS647425	SIGVARIS Cotton - Knee - 18-21mmHg XL Norm Nature Closed Toe	Natural
PTS647420	SIGVARIS Cotton - Knee - 18-21mmHg XL Norm Black Open Toe	Black
PTS647415	SIGVARIS Cotton - Knee - 18-21mmHg XL Norm Black Closed Toe	Black
PTS647410	SIGVARIS Cotton - Knee - 18-21mmHg XL Long Nature Open Toe	Natural
PTS647405	SIGVARIS Cotton - Knee - 18-21mmHg XL Long Nature Closed Toe	Natural
PTS647400	SIGVARIS Cotton - Knee - 18-21mmHg XL Long Black Open Toe	Black
PTS647335	SIGVARIS Cotton - Knee - 18-21mmHg M Norm Black Closed Toe	Black
PTS647390	SIGVARIS Cotton - Knee - 18-21mmHg S Norm Nature Open Toe	Natural
PTS647340	SIGVARIS Cotton - Knee - 18-21mmHg M Norm Black Open Toe	Black
PTS647380	SIGVARIS Cotton - Knee - 18-21mmHg S Norm Black Open Toe	Black
PTS647375	SIGVARIS Cotton - Knee - 18-21mmHg S Norm Black Closed Toe	Black
PTS647370	SIGVARIS Cotton - Knee - 18-21mmHg S Long Nature Open Toe	Natural
PTS647365	SIGVARIS Cotton - Knee - 18-21mmHg S Long Nature Closed Toe	Natural
PTS647360	SIGVARIS Cotton - Knee - 18-21mmHg S Long Black Open Toe	Black
PTS647355	SIGVARIS Cotton - Knee - 18-21mmHg S Long Black Closed Toe	Black
PTS647350	SIGVARIS Cotton - Knee - 18-21mmHg M Norm Nature Open Toe	Natural
PTS647230	SIGVARIS Cotton - Knee - 18-21mmHg XS PLUS Norm Nature Open Toe	Natural

Code	Display Name	Colour
PTS647395	SIGVARIS Cotton - Knee - 18-21mmHg XL Long Black Closed Toe	Black
PTS647075	SIGVARIS Cotton - Knee - 18-21mmHg M PLUS Long Black Closed Toe	Black
PTS647240	SIGVARIS Cotton - Knee - 18-21mmHg XXS PLUS Long Black Open Toe	Black
PTS647120	SIGVARIS Cotton - Knee - 18-21mmHg S PLUS Long Black Open Toe	Black
PTS647115	SIGVARIS Cotton - Knee - 18-21mmHg S PLUS Long Black Closed Toe	Black
PTS647110	SIGVARIS Cotton - Knee - 18-21mmHg M PLUS Norm Nature Open Toe	Natural
PTS647105	SIGVARIS Cotton - Knee - 18-21mmHg M PLUS Norm Nature Closed Toe	Natural
PTS647100	SIGVARIS Cotton - Knee - 18-21mmHg M PLUS Norm Black Open Toe	Black
PTS647095	SIGVARIS Cotton - Knee - 18-21mmHg M PLUS Norm Black Closed Toe	Black
PTS647090	SIGVARIS Cotton - Knee - 18-21mmHg M PLUS Long Nature Open Toe	Natural
PTS647130	SIGVARIS Cotton - Knee - 18-21mmHg S PLUS Long Nature Open Toe	Natural
PTS647080	SIGVARIS Cotton - Knee - 18-21mmHg M PLUS Long Black Open Toe	Black
PTS647135	SIGVARIS Cotton - Knee - 18-21mmHg S PLUS Norm Black Closed Toe	Black
PTS647070	SIGVARIS Cotton - Knee - 18-21mmHg L PLUS Norm Nature Open Toe	Natural
PTS647065	SIGVARIS Cotton - Knee - 18-21mmHg L PLUS Norm Nature Closed Toe	Natural
PTS647060	SIGVARIS Cotton - Knee - 18-21mmHg L PLUS Norm Black Open Toe	Black
PTS647055	SIGVARIS Cotton - Knee - 18-21mmHg L PLUS Norm Black Closed Toe	Black
PTS647050	SIGVARIS Cotton - Knee - 18-21mmHg L PLUS Long Nature Open Toe	Natural

Code	Display Name	Colour
PTS647045	SIGVARIS Cotton - Knee - 18-21mmHg L PLUS Long Nature Closed Toe	Natural
PTS647040	SIGVARIS Cotton - Knee - 18-21mmHg L PLUS Long Black Open Toe	Black
PTS647035	SIGVARIS Cotton - Knee - 18-21mmHg L PLUS Long Black Closed Toe	Black
PTS647085	SIGVARIS Cotton - Knee - 18-21mmHg M PLUS Long Nature Closed Toe	Natural
PTS647180	SIGVARIS Cotton - Knee - 18-21mmHg XL PLUS Norm Black Open Toe	Black
PTS647445	SIGVARIS Cotton - Knee - 18-21mmHg XS Long Nature Closed Toe	Natural
PTS647225	SIGVARIS Cotton - Knee - 18-21mmHg XS PLUS Norm Nature Closed Toe	Natural
PTS647220	SIGVARIS Cotton - Knee - 18-21mmHg XS PLUS Norm Black Open Toe	Black
PTS647215	SIGVARIS Cotton - Knee - 18-21mmHg XS PLUS Norm Black Closed Toe	Black
PTS647210	SIGVARIS Cotton - Knee - 18-21mmHg XS PLUS Long Nature Open Toe	Natural
PTS647205	SIGVARIS Cotton - Knee - 18-21mmHg XS PLUS Long Nature Closed Toe	Natural
PTS647200	SIGVARIS Cotton - Knee - 18-21mmHg XS PLUS Long Black Open Toe	Black
PTS647195	SIGVARIS Cotton - Knee - 18-21mmHg XS PLUS Long Black Closed Toe	Black
PTS647125	SIGVARIS Cotton - Knee - 18-21mmHg S PLUS Long Nature Closed Toe	Natural
PTS647185	SIGVARIS Cotton - Knee - 18-21mmHg XL PLUS Norm Nature Closed Toe	Natural
PTS647235	SIGVARIS Cotton - Knee - 18-21mmHg XXS PLUS Long Black Closed Toe	Black
PTS647175	SIGVARIS Cotton - Knee - 18-21mmHg XL PLUS Norm Black Closed Toe	Black
PTS647170	SIGVARIS Cotton - Knee - 18-21mmHg XL PLUS Long Nature Open Toe	Natural

Code	Display Name	Colour
PTS647165	SIGVARIS Cotton - Knee - 18-21mmHg XL PLUS Long Nature Closed Toe	Natural
PTS647160	SIGVARIS Cotton - Knee - 18-21mmHg XL PLUS Long Black Open Toe	Black
PTS647155	SIGVARIS Cotton - Knee - 18-21mmHg XL PLUS Long Black Closed Toe	Black
PTS647150	SIGVARIS Cotton - Knee - 18-21mmHg S PLUS Norm Nature Open Toe	Natural
PTS647145	SIGVARIS Cotton - Knee - 18-21mmHg S PLUS Norm Nature Closed Toe	Natural
PTS647140	SIGVARIS Cotton - Knee - 18-21mmHg S PLUS Norm Black Open Toe	Black
PTS647190	SIGVARIS Cotton - Knee - 18-21mmHg XL PLUS Norm Nature Open Toe	Natural
PTS647700	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Norm Black Open Toe	Black
PTS647745	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Long Nature Closed Toe	Natural
PTS647740	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Long Black Open Toe	Black
PTS647735	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Long Black Closed Toe	Black
PTS647730	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Short Nature Open Toe	Natural
PTS647725	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Short Nature Closed Toe	Natural
PTS647720	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Short Black Open Toe	Black
PTS647715	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Short Black Closed Toe	Black
PTS647435	SIGVARIS Cotton - Knee - 18-21mmHg XS Long Black Closed Toe	Black
PTS647705	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Norm Nature Closed Toe	Natural
PTS647760	SIGVARIS Cotton - Knee - 23-32mmHg S Long Black Open Toe	Black

Code	Display Name	Colour
PTS647695	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Norm Black Closed Toe	Black
PTS647690	SIGVARIS Cotton - Knee - 23-32mmHg M Long Nature Open Toe	Natural
PTS647685	SIGVARIS Cotton - Knee - 23-32mmHg M Long Nature Closed Toe	Natural
PTS647680	SIGVARIS Cotton - Knee - 23-32mmHg M Long Black Open Toe	Black
PTS647675	SIGVARIS Cotton - Knee - 23-32mmHg M Long Black Closed Toe	Black
PTS647670	SIGVARIS Cotton - Knee - 23-32mmHg M Short Nature Open Toe	Natural
PTS647665	SIGVARIS Cotton - Knee - 23-32mmHg M Short Nature Closed Toe	Natural
PTS647660	SIGVARIS Cotton - Knee - 23-32mmHg M Short Black Open Toe	Black
PTS647710	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Norm Nature Open Toe	Natural
PTS647800	SIGVARIS Cotton - Knee - 23-32mmHg S Short Black Open Toe	Black
PTS648685	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg M Long Nature Closed Toe	Natural
PTS647845	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Norm Nature Closed Toe	Natural
PTS647840	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Norm Black Open Toe	Black
PTS647835	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Norm Black Closed Toe	Black
PTS647830	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Long Nature Open Toe	Natural
PTS647825	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Long Nature Closed Toe	Natural
PTS647820	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Long Black Open Toe	Black
PTS647815	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Long Black Closed Toe	Black

Code	Display Name	Colour
PTS647750	SIGVARIS Cotton - Knee - 23-32mmHg M PLUS Long Nature Open Toe	Natural
PTS647805	SIGVARIS Cotton - Knee - 23-32mmHg S Short Nature Closed Toe	Natural
PTS647755	SIGVARIS Cotton - Knee - 23-32mmHg S Long Black Closed Toe	Black
PTS647795	SIGVARIS Cotton - Knee - 23-32mmHg S Short Black Closed Toe	Black
PTS647790	SIGVARIS Cotton - Knee - 23-32mmHg S Norm Nature Open Toe	Natural
PTS647785	SIGVARIS Cotton - Knee - 23-32mmHg S Norm Nature Closed Toe	Natural
PTS647780	SIGVARIS Cotton - Knee - 23-32mmHg S Norm Black Open Toe	Black
PTS647775	SIGVARIS Cotton - Knee - 23-32mmHg S Norm Black Closed Toe	Black
PTS647770	SIGVARIS Cotton - Knee - 23-32mmHg S Long Nature Open Toe	Natural
PTS647765	SIGVARIS Cotton - Knee - 23-32mmHg S Long Nature Closed Toe	Natural
PTS647645	SIGVARIS Cotton - Knee - 23-32mmHg M Norm Nature Closed Toe	Natural
PTS647810	SIGVARIS Cotton - Knee - 23-32mmHg S Short Nature Open Toe	Natural
PTS647490	SIGVARIS Cotton - Knee - 18-21mmHg XXS Long Nature Open Toe	Natural
PTS647655	SIGVARIS Cotton - Knee - 23-32mmHg M Short Black Closed Toe	Black
PTS647535	SIGVARIS Cotton - Knee - 23-32mmHg L Norm Black Closed Toe	Black
PTS647530	SIGVARIS Cotton - Knee - 23-32mmHg L Long Nature Open Toe	Natural
PTS647525	SIGVARIS Cotton - Knee - 23-32mmHg L Long Nature Closed Toe	Natural
PTS647520	SIGVARIS Cotton - Knee - 23-32mmHg L Long Black Open Toe	Black

Code	Display Name	Colour
PTS647515	SIGVARIS Cotton - Knee - 23-32mmHg L Long Black Closed Toe	Black
PTS647510	SIGVARIS Cotton - Knee - 18-21mmHg XXS Norm Nature Open Toe	Natural
PTS647505	SIGVARIS Cotton - Knee - 18-21mmHg XXS Norm Nature Closed Toe	Natural
PTS647545	SIGVARIS Cotton - Knee - 23-32mmHg L Norm Nature Closed Toe	Natural
PTS647495	SIGVARIS Cotton - Knee - 18-21mmHg XXS Norm Black Closed Toe	Black
PTS647550	SIGVARIS Cotton - Knee - 23-32mmHg L Norm Nature Open Toe	Natural
PTS647485	SIGVARIS Cotton - Knee - 18-21mmHg XXS Long Nature Closed Toe	Natural
PTS647480	SIGVARIS Cotton - Knee - 18-21mmHg XXS Long Black Open Toe	Black
PTS647475	SIGVARIS Cotton - Knee - 18-21mmHg XXS Long Black Closed Toe	Black
PTS647470	SIGVARIS Cotton - Knee - 18-21mmHg XS Norm Nature Open Toe	Natural
PTS647465	SIGVARIS Cotton - Knee - 18-21mmHg XS Norm Nature Closed Toe	Natural
PTS647460	SIGVARIS Cotton - Knee - 18-21mmHg XS Norm Black Open Toe	Black
PTS647455	SIGVARIS Cotton - Knee - 18-21mmHg XS Norm Black Closed Toe	Black
PTS647450	SIGVARIS Cotton - Knee - 18-21mmHg XS Long Nature Open Toe	Natural
PTS647500	SIGVARIS Cotton - Knee - 18-21mmHg XXS Norm Black Open Toe	Black
PTS647595	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Short Black Closed Toe	Black
PTS647855	SIGVARIS Cotton - Knee - 23-32mmHg S PLUS Short Black Closed Toe	Black
PTS647640	SIGVARIS Cotton - Knee - 23-32mmHg M Norm Black Open Toe	Black

Code	Display Name	Colour
PTS647635	SIGVARIS Cotton - Knee - 23-32mmHg M Norm Black Closed Toe	Black
PTS647630	SIGVARIS Cotton - Knee - 23-32mmHg L Short Nature Open Toe	Natural
PTS647625	SIGVARIS Cotton - Knee - 23-32mmHg L Short Nature Closed Toe	Natural
PTS647620	SIGVARIS Cotton - Knee - 23-32mmHg L Short Black Open Toe	Black
PTS647615	SIGVARIS Cotton - Knee - 23-32mmHg L Short Black Closed Toe	Black
PTS647610	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Short Nature Open Toe	Natural
PTS647540	SIGVARIS Cotton - Knee - 23-32mmHg L Norm Black Open Toe	Black
PTS647600	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Short Black Open Toe	Black
PTS647650	SIGVARIS Cotton - Knee - 23-32mmHg M Norm Nature Open Toe	Natural
PTS647590	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Norm Nature Open Toe	Natural
PTS647585	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Norm Nature Closed Toe	Natural
PTS647580	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Norm Black Open Toe	Black
PTS647575	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Norm Black Closed Toe	Black
PTS647570	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Long Nature Open Toe	Natural
PTS647565	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Long Nature Closed Toe	Natural
PTS647560	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Long Black Open Toe	Black
PTS647555	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Long Black Closed Toe	Black
PTS647605	SIGVARIS Cotton - Knee - 23-32mmHg L PLUS Short Nature Closed Toe	Natural

Code	Display Name	Colour
PTS649005	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XXS Short Nature Closed Toe	Natural
PTS649065	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg M Short Nature Closed Toe	Natural
PTS649060	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg M Norm Nature Open Toe	Natural
PTS649055	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg M Norm Nature Closed Toe	Natural
PTS649050	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg M Long Nature Open Toe	Natural
PTS649045	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg M Long Nature Closed Toe	Natural
PTS649040	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg L Short Nature Open Toe	Natural
PTS649035	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg L Short Nature Closed Toe	Natural
PTS649030	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg L Norm Nature Open Toe	Natural
PTS649025	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg L Norm Nature Closed Toe	Natural
PTS649020	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg L Long Nature Open Toe	Natural
PTS648940	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XL Norm Nature Open Toe	Natural
PTS649010	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XXS Short Nature Open Toe	Natural
PTS649080	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg S Long Nature Open Toe	Natural
PTS649000	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XXS Norm Nature Open Toe	Natural
PTS648995	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XXS Norm Nature Closed Toe	Natural
PTS648990	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XXS Long Nature Open Toe	Natural
PTS648980	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XS Short Nature Open Toe	Natural

Code	Display Name	Colour
PTS648970	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XS Norm Nature Open Toe	Natural
PTS648965	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XS Norm Nature Closed Toe	Natural
PTS648960	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XS Long Nature Open Toe	Natural
PTS648955	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XS Long Nature Closed Toe	Natural
PTS648950	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XL Short Nature Open Toe	Natural
PTS648945	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XL Short Nature Closed Toe	Natural
PTS649015	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg L Long Nature Closed Toe	Natural
PTS649135	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XS Long Nature Closed Toe	Natural
PTS649190	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XXS Short Nature Open Toe	Natural
PTS649185	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XXS Short Nature Closed Toe	Natural
PTS649180	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XXS Norm Nature Open Toe	Natural
PTS649175	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XXS Norm Nature Closed Toe	Natural
PTS649170	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XXS Long Nature Open Toe	Natural
PTS649165	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XXS Long Nature Closed Toe	Natural
PTS649160	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XS Short Nature Open Toe	Natural
PTS649155	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XS Short Nature Closed Toe	Natural
PTS649150	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XS Norm Nature Open Toe	Natural
PTS649070	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg M Short Nature Open Toe	Natural

Code	Display Name	Colour
PTS649140	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XS Long Nature Open Toe	Natural
PTS649075	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg S Long Nature Closed Toe	Natural
PTS649130	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XL Short Nature Open Toe	Natural
PTS649125	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XL Short Nature Closed Toe	Natural
PTS649120	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XL Norm Nature Open Toe	Natural
PTS649115	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XL Norm Nature Closed Toe	Natural
PTS649110	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XL Long Nature Open Toe	Natural
PTS649105	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XL Long Nature Closed Toe	Natural
PTS649100	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg S Short Nature Open Toe	Natural
PTS649095	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg S Short Nature Closed Toe	Natural
PTS649090	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg S Norm Nature Open Toe	Natural
PTS649085	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg S Norm Nature Closed Toe	Natural
PTS648975	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XS Short Nature Closed Toe	Natural
PTS649145	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 23-32mmHg XS Norm Nature Closed Toe	Natural
PTS648740	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg S Short Nature Open Toe	Natural
PTS648805	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XXS Long Nature Closed Toe	Natural
PTS648800	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XS Short Nature Open Toe	Natural
PTS648795	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XS Short Nature Closed Toe	Natural

Code	Display Name	Colour
PTS648790	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XS Norm Nature Open Toe	Natural
PTS648785	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XS Norm Nature Closed Toe	Natural
PTS648780	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XS Long Nature Open Toe	Natural
PTS648935	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XL Norm Nature Closed Toe	Natural
PTS648770	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XL Short Nature Open Toe	Natural
PTS648985	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XXS Long Nature Closed Toe	Natural
PTS648760	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XL Norm Nature Open Toe	Natural
PTS648755	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XL Norm Nature Closed Toe	Natural
PTS648810	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XXS Long Nature Open Toe	Natural
PTS648745	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XL Long Nature Closed Toe	Natural
PTS648775	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XS Long Nature Closed Toe	Natural
PTS648735	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg S Short Nature Closed Toe	Natural
PTS648730	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg S Norm Nature Open Toe	Natural
PTS648725	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg S Norm Nature Closed Toe	Natural
PTS648720	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg S Long Nature Open Toe	Natural
PTS648715	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg S Long Nature Closed Toe	Natural
PTS648710	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg M Short Nature Open Toe	Natural
PTS648705	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg M Short Nature Closed Toe	Natural

Code	Display Name	Colour
PTS648700	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg M Norm Nature Open Toe	Natural
PTS648695	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg M Norm Nature Closed Toe	Natural
PTS648675	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg L Short Nature Closed Toe	Natural
PTS648690	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg M Long Nature Open Toe	Natural
PTS648750	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XL Long Nature Open Toe	Natural
PTS648905	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg S Norm Nature Closed Toe	Natural
PTS648930	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XL Long Nature Open Toe	Natural
PTS648925	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg XL Long Nature Closed Toe	Natural
PTS648920	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg S Short Nature Open Toe	Natural
PTS648765	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XL Short Nature Closed Toe	Natural
PTS648910	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg S Norm Nature Open Toe	Natural
PTS648815	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XXS Norm Nature Closed Toe	Natural
PTS648900	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg S Long Nature Open Toe	Natural
PTS648895	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg S Long Nature Closed Toe	Natural
PTS648890	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg M Short Nature Open Toe	Natural
PTS648885	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg M Short Nature Closed Toe	Natural
PTS648880	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg M Norm Nature Open Toe	Natural
PTS648835	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg L Long Nature Closed Toe	Natural

Code	Display Name	Colour
PTS648915	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg S Short Nature Closed Toe	Natural
PTS648875	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg M Norm Nature Closed Toe	Natural
PTS648820	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XXS Norm Nature Open Toe	Natural
PTS648830	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XXS Short Nature Open Toe	Natural
PTS648840	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg L Long Nature Open Toe	Natural
PTS648845	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg L Norm Nature Closed Toe	Natural
PTS648850	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg L Norm Nature Open Toe	Natural
PTS648855	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg L Short Nature Closed Toe	Natural
PTS648860	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg L Short Nature Open Toe	Natural
PTS648865	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg M Long Nature Closed Toe	Natural
PTS648870	SIGVARIS Cotton - Thigh w/Grip Top PLUS - 18-21mmHg M Long Nature Open Toe	Natural
PTS648825	SIGVARIS Cotton - Thigh w/Grip Top - 23-32mmHg XXS Short Nature Closed Toe	Natural