

Jobst Bella Lite - Combination Sleeve/Glove (20 - 30 mmHg)



- Soft and air-permeable
- Lightweight knit construction
- Fits a wide range of hand and arm shapes
- Fits both the left or right arm/hand
- The combination sleeve/glove can be worn separately together

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information



Code	Display Name	Elbow Circumference	Wrist Circumference
FBS303484	Jobst Bella Lite - Ready To Wear - Combined 20-30mmHg Sml/Reg	220 - 400 mm	180 - 230 mm
FBS303485	Jobst Bella Lite - Ready To Wear - Combined 20-30mmHg Med/Reg	270 -450 mm	230 - 290 mm
FBS303486	Jobst Bella Lite - Ready To Wear - Combined 20-30mmHg Lge/Reg	330 - 510 mm	290 - 340 mm
FBS303487	Jobst Bella Lite - Ready To Wear - Combined 20-30mmHg Sml/Long	220 - 400 mm	180 - 230 mm
FBS303488	Jobst Bella Lite - Ready To Wear - Combined 20-30mmHg Med/Long	270 -450 mm	230 - 290 mm
FBS303489	Jobst Bella Lite - Ready To Wear - Combined 20-30mmHg Lge/Long	330 - 510 mm	290 - 340 mm