

Aspire Montana Single Action Maxi Lift Recline Chair



- The Montana Lift Recline Chair provides a large seating surface, offering improved comfort for users who need that little bit of extra room.
- Single motor lift recline chair
- Designed for larger users
- Head & armrest covers for improved longevity
- Superior wipe clean microsuede upholstery
- Part of the Aspire BRONZE range of lift recline chairs
- Available in Bronze Microsuede Leather-look
- Looking to protect your Electric Lift Recline Chair? Try our [Universal Chair Protector](#).

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_180903180_0103_1_2

Click [here](#) for more information

Compliance

ISO13934-1 / ISO13937-1 / BS6807
/ BS7175

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load
CHP227400	Chair Powerlift Recline - Aspire Montana Bariatric Recline - Bronze - Single Motor - 220kg	2659	1	530 mm	630 mm	660 mm	510 mm	220 Kg