

DonJoy Universal Arm Sling



- Durable cotton/poly material
- Ideal for immobilisation and support of shoulder and elbow joints
- Contact closure adjustment on envelope accommodates sizing for most patients
- Double O rings and contact closure shoulder strap for proper adjustment
- Available with or without fully padded strap

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Length
FBS281990	DonJoy Universal Arm Sling	457 mm