

Aidacare Anaesthesia Cart



- The Aidacare Anaesthesia Cart is the flexible solution to your workflow needs
- Aidacare carts are designed to offer an attractive, practical and durable cart with multiple configurations to enable the efficient organisation and delivery of medication and equipment
- Robust panels to minimise dents and scratches
- Combined with its sleek lines, strong construction and easy manoeuvrability, it offers a fully customisable solution to meet individual needs
- Easy roll 125 mm castors
- Packed with standard accessories, it allows for choice of cart height, draw depth, locking option and storage configuration to perfectly compliment your workflow needs
- **Standard Inclusions:** Accessory Tower, Drawer Dividers, Extended Sliding Surface and Waste Bin with Lid
- Choose from an array of accessories to optimise the cart for your facility's treatment needs
- Lock System Options:
 - **Break-away lock with Plastic Seal:** The simple breakaway lock offers easy access, combined with essential security
 - **Digital lock:** A user friendly system where the cart is locked and unlocked using a pre-set code ensuring ease of access and the security required in clinical environments. The digital lock is robust and waterproof
- Standard drawer configuration detailed below
- **Small Anaesthesia Cart:**
 - 1 small 3' drawer (7 cm)
 - 2 medium 6' drawers (16 cm)
 - 1 large 9' drawer (25 cm)
- **Medium Anaesthesia Cart:**



Accessories

Click on the cart icon to display the range of optional accessories

- 2 small 3' drawers (7 cm)
- 2 medium 6' drawers (16 cm)
- 1 large 9' drawer (25 cm)
- **Large Anaesthesia Cart:**
 - 1 small 3' drawer (7 cm)
 - 3 medium 6' drawers (16 cm)
 - 1 large 9' drawer (25cm)

Warranty Information

- Structural: 5 Years Warranty
- Digital Lock: 2 Years Warranty
- Accessories: 1 Year Warranty

Code	Display Name	Total Length	Total Width	Total Height	Colour
WAT909660	Aidacare Anaesthesia Cart - Small - Break-away Lock	790 mm	610 mm	900 mm	Green
WAT909630	Aidacare Anaesthesia Cart - Small - Digital Lock	790 mm	610 mm	900 mm	Green
WAT909670	Aidacare Anaesthesia Cart - Medium - Break-away Lock	790 mm	610 mm	1000 mm	Green
WAT909640	Aidacare Anaesthesia Cart - Medium - Digital Lock	790 mm	610 mm	1000 mm	Green
WAT909680	Aidacare Anaesthesia Cart - Large - Break-away Lock	790 mm	610 mm	1100 mm	Green
WAT909650	Aidacare Anaesthesia Cart - Large - Digital Lock	790 mm	610 mm	1100 mm	Green