

## Acute Bedside Cabinet



- The Acute Bedside Cabinet offers versatile storage with simple and functional design.
- Features two secure lockable drawers to keep medication and valuables secure.
- Purpose-built design to meet the rigours of hospital and aged care applications.
- The slimline design will fit into a facility, bedroom or common area easily.
- Rear towel rail offers a practical and tidy solution for compact living.
- Additional drawers provide ample space for bedside essentials and personal items
- The open top provides space for a lamp, book and smart devices.
- Select from our range of laminate colour options.



Code	Display Name	Total Length	Total Width	Total Height
BEA826380	Acute Bedside Cabinet - Laminate	525 mm	465 mm	840 mm