

Rosie Cool Highback Chair



- Sleek and modern armchair that suits a contemporary interior
- An eye-catching design with smooth, straight lines
- Sturdy, comfortable back and cushioned seat offer a welcoming place to sit and relax
- Narrow solid timber legs with a walnut finish provide a solid base
- Quality construction ensures this armchair will withstand daily use
- Fabric colour/material to be advised at time of order



Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
FUL052200	Rosie Cool Highback Chair	760mm	660mm	1110mm