

Bath Transfer Bench



- Provides easy and safe access as well as a stable seating area for bathing and showering.
- Adjustable height aluminium frame to suit individual user requirements.
- Plastic moulded textured top features holes to drain water away from the seating area.
- This transfer bench has a sliding plastic seat and backrest to assist users who have difficulty with manually moving across a bench.
- A velcro lap belt securely holds the user, and the backrest is useful for those with reduced sitting balance.

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
BTS128600	Bath Transfer Bench with Sliding Seat	100 kg	340 mm	520 mm	460 - 560 mm	460 - 560 mm