

## Bath Transfer Bench



- Provides easy and safe access as well as a stable seating area for bathing and showering.
- Adjustable height aluminium frame to suit individual user requirements.
- Plastic moulded textured top features holes to drain water away from the seating area.
- This transfer bench has a sliding plastic seat and backrest to assist users who have difficulty with manually moving across a bench.
- A velcro lap belt securely holds the user, and the backrest is useful for those with reduced sitting balance.

- ✓ Contracted DVA product
- ✓ Contract approved

Code	Display Name	RRP	Box Quantity	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
BTS128600	Bath Transfer Bench with Sliding Seat	679	1	100 kg	340 mm	520 mm	460 - 560 mm	460 - 560 mm