

## **Panamatic Midi**





- The Panamatic Midi Flusher Thermal Sanitiser is the cold water cousin of the
  Optima 2 and is designed for low-medium workload situations
- Improved Wash Performance (Cold water) An increased number of more efficient cleaning nozzles and optimum water jet pressure result in performance levels that make the Midi outstanding for its size
- Reduced Cycle-Time Better washing and temperature control has enabled improved overall cycle-time for quicker turnaround
- More Efficient Water and Power With improved performance and reduced cycle time, both power and water consumption have been reduced. Cold water flush / wash - steam cycle - cool-down rinse
- Hands-free Opening and Cycle Start A foot operated opening mechanism and hands-free close and cycle start means direct contact with the machine has been reduced to a minimum, helping to reduce potential cross infection
- MicrobeSafe Antimicrobial Lid DDC™ unique lid incorporates Silver
  Nanotechnology to help stop bacterial growth and reduce cross infection
- Large Chamber Capacity and Easier Loading The wash chamber and larger opening allow easier loading of a wide range of utensils. The chamber can process 2 bedpans and 2 urine bottles per load or various combinations
- Electrical supply: Single Phase 230 V, 13 AMP
- Electrical output/input: Motor 0.425 KW, Heater 2.5 KW
- Water supply: 15 mm Cold water 0.5 10 bar
- Soil connection: 110 mm connections supplied
- Operation cycle time: Approx 5 mins / 30 ml EcoWash
- Steam disinfection: Exceeds Ao 700, 85° x 75 secs

## **Warranty Information**

• 1 Year



## **Compliance**

ISO 15883/AS 4187:2019



Code Display Name

UPM100310 DDC Bedpan Washer Sanitiser - Midi