

## Permobil Wheelchair Elastic Chest Strap



- For anytime positioning at home, in clinic, fittings, motor vehicles, wheeling over rough streets, or for sports, the Universal Elastic Strap is the stable solution.
- Strong elastic strap can be used on the chest, abdomen, lap or legs, and has an easy-grasp end tab with a thumbhole for users with limited hand function – a versatile support that’s more comfortable than traditional binders or corsets.
- Can be wrapped around a wheelchair back, tubing, motor vehicle seat or recliner to ensure that you get the support you need where you need it.



| Code      | Display Name   |
|-----------|--|
| WCA860539 | Wheelchair Chest Strap - Elastic - Small - 142cm x 8cm |
| WCA860540 | Wheelchair Chest Strap - Elastic - Lge - 142cm x 11cm  |
| WCA860541 | Wheelchair Chest Strap - Elastic - XL - 142cm x 15cm   |