

Permobil Wheelchair Elastic Chest Strap



- For anytime positioning at home, in clinic, fittings, motor vehicles, wheeling over rough streets, or for sports, the Universal Elastic Strap is the stable solution.
- Strong elastic strap can be used on the chest, abdomen, lap or legs, and has an easy-grasp end tab with a thumbhole for users with limited hand function – a versatile support that’s more comfortable than traditional binders or corsets.
- Can be wrapped around a wheelchair back, tubing, motor vehicle seat or recliner to ensure that you get the support you need where you need it.

 Contract approved

Support @ Home

Function: Wheelchair accessories

Category: Prescribed

Click [here](#) for more information

Code	Display Name
WCA860539	Wheelchair Chest Strap - Elastic - Small - 142cm x 8cm
WCA860540	Wheelchair Chest Strap - Elastic - Lge - 142cm x 11cm
WCA860541	Wheelchair Chest Strap - Elastic - XL - 142cm x 15cm