

Forearm Crutches - Custom



- Reinforced aluminium construction designed for users with long term needs
- Absence of adjustment holes makes the crutch 'silent' in use
- Supplied in one standard length and may be cut to the required length
- Please specify handle to floor height at time of order
- Ferrule: 25 mm
- Product Weight: 800 g (each)

Warranty Information

- 1 Year

 Contract approved

Support @ Home

Function: Walking, manipulated by one arm

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_120309088_0105_1_2

Walking Aids -

Composite:05_416_0105_1_2

Click [here](#) for more information

Compliance

ISO 11199-2:2005

Code	Display Name	Safe Working Load	Unit Weight
WAC688400	Forearm Crutches - Custom	190 kg	800 g (each)