

MEDI K.400 EPICOMED ELASTIC ELBOW



- The compressive knitted fabric and silicone insert exert a massage effect which promotes the healing process
- Specific areas of pressure to reduce pain and relieve tension at the tendon insert Application:
- Lateral and medial epicondylitis (tennis elbow, golfer's elbow)
- Joint effusion and swelling due to arthritis and osteoarthritis
- Chronic post-operative and post-traumatic inflammation

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS307280	MEDI K.400 EPICOMED ELASTIC ELBOW - SIZE I
FBS307290	MEDI K.400 EPICOMED ELASTIC ELBOW - SIZE II
FBS307300	MEDI K.400 EPICOMED ELASTIC ELBOW - SIZE III
FBS307310	MEDI K.400 EPICOMED ELASTIC ELBOW - SIZE IV
FBS307320	MEDI K.400 EPICOMED ELASTIC ELBOW - SIZE V
FBS307330	MEDI K.400 EPICOMED ELASTIC ELBOW - SIZE VI