

Digi Flex Exercise System



- Develops isolated finger strength, flexibility and coordination
- Ergonomic base is easy to grip and made with non-latex material
- Available resistances range from .34 kg to 5.9kg (colour-coded)

✔ Contract approved

Support @ Home

Function: Not on S@H

Category: Not on S@H

Click [here](#) for more information

Code	Display Name
FPP298700	Digi Flex Exercise System - Yellow 0.7kg
FPP298800	Digi Flex Exercise System - Red 1.4kg
FPP298900	Digi Flex Exercise System - Green 2.3kg
FPP299100	Digi Flex Exercise System - Black 4.1kg