

Padded Transfer Belt



- Assists the carer in patient mobility and transport, padded for user comfort
- Reduces risk of strain or injury to carer
- Adjustable and with nylon buckle to suit user requirements

Warranty Information

- 1 Year

- ✓ Contracted DVA product
- ✓ Contract approved

ndis

Item Code:

05_123115159_0105_1_2

Minor Transfer –

Composite:05_430_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	RRP	Box Quantity	Total Length
TSE673200	Padded Transfer Belt - Small	253	1	750 - 850 mm
TSE673190	Padded Transfer Belt - Medium	253	1	850 - 950 mm
TSE673180	Padded Transfer Belt - Large	253	1	950 - 1050 mm