

Padded Transfer Belt



- Assists the carer in patient mobility and transport, padded for user comfort
- Reduces risk of strain or injury to carer
- Adjustable and with nylon buckle to suit user requirements

Warranty Information

- 1 Year

- ✔ Contracted DVA product
- ✔ Contract approved

Support @ Home

Function: Changing body position

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123115159_0105_1_2

Minor Transfer –

Composite:05_430_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length	Material
TSE673200	Padded Transfer Belt - Small	750 - 850 mm	Fabric
TSE673190	Padded Transfer Belt - Medium	850 - 950 mm	Fabric

Code	Display Name	Total Length	Material
TSE673180	Padded Transfer Belt - Large	950 - 1050 mm	Fabric