

Derik Bedside Cabinet



- The Derik Bedside Cabinet allows easy access to bedtime essentials with its clean lines and simple aesthetic
- The simple silhouette features two drawers and a underdraw compartment for ample storage
- Raised safety edging on the cabinet top prevents items from sliding off
- Rear towel rail offers a practical and tidy solution for compact living
- Made with quality materials to ensure lasting functionality
- Slimline handles and smooth drawer action offer function and practicality
- The open top provides space for a lamp, books and smart devices
- Optional hutch provides added shelf space along with practical storage for books and magazines
- Select from our range of laminate colour options



| Code | Display Name | Total Length | Total Width | Total Height |
|-----------|----------------------------------|--------------|-------------|--------------|
| BEA826340 | Derik Bedside Cabinet - Laminate | 450 mm | 465 mm | 810 mm |