

Active Air 8 Alternating Mattress Replacement



The Active Air 8 is a dynamic mattress replacement suitable for users at “High Risk” of Pressure Injury. Featuring AUTODETECT, and AUTOPROTECT, together with CELL-IN-CELL technology, the Active Air 8 provides the highest levels of care into Acute and Long Term care environments. The Active Air 8’s innovative features ensure lowest interface pressures with an easy-to-use digital pump.

- Digital Pump – The Easy-to-read LCD screen provides clear and precise information on cycle time, operating mode and fault codes to assist with troubleshooting
- AutoLock – Control panel automatically locks after 30 seconds, preventing accidental mode changes
- [Watch instructional video to setting up your mattress](#)
- AutoDetect – Based on the user’s weight and size, appropriate pressure output is automatically set, ensuring optimal care and protection
- AutoProtect – When the head of the bed is raised to 45° or more, air pressure automatically increases to prevent bottoming out
- Continuous Low Pressure (CLP) – When activated, this mode increases user immersion and envelopment
- 24hr Transport Mode – In the event of a power failure, this feature provides a security buffer allowing patients to remain on an inflated surface until power is restored
- Sleep Quality - Innovative pump technology reduces noise and vibration, for a quieter sleep environment
- Infection Control - Welded seams and zip protection prevent fluid and bacterial ingress into the system
- Skin Microclimate - The multi-stretch cover reduces shear and offers a high Moisture Vapour transfer rate (MVTR)
- Workplace Health & Safety - The power cord is retained within the overlay

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_043306005_0103_1_2

Power Bed – Composite:

05_423_0103_1_2

Click [here](#) for more information

Compliance

IEC 60601-1-2

Accessories

Click on the cart icon to display the range of optional accessories

eliminating its potential as a trip hazard

- Pressure Care Outcomes - Intelligent cell design maximises the duration of pressure reduction in the alternation cycle
- Laundry & Hygiene - Machine washable cover at 95° to achieve superior sanitising

Warranty Information

- Pump: 2 Years
- Replacement: 2 Years

| Code | Display Name | RRP | Total Length | Total Width | Total Height | Safe Working Load | Pressure Care Rating |
|-----------|---|------|--------------|-------------|--------------|-------------------|----------------------|
| PTM639300 | Active Air 8 Alternating Mattress Replacement - Single | 3539 | 2000 mm | 900 mm | 200 mm | 30 - 200 kg | High Risk |
| PTM639310 | Active Air 8 Alternating Mattress Replacement - King Single | 4919 | 2000 mm | 1050 mm | 200 mm | 30 - 200 kg | High Risk |
| PTM639000 | Rapid Inflation Turbo Pump | 432 | - | - | - | - | - |
| PTM639020 | Carry Bag - Active Air 8/8KS/BARI 10KS Mattress | 103 | - | - | - | - | - |