

Affinity Grab Handle / Bed Rail



The Affinity Grab Handle/Bed Rail provides reliable and stable support for users during bed-related activities, including sitting, standing, and transitioning between lying and sitting positions. It anchors securely to the underlying bed frame and removes the potential for shifting during a transfer.



- Support frames slide beneath the mattress and attach securely to the bed base, providing additional stability during use
- Safety straps fasten beneath the bed to prevent shifting and eliminate trip hazards
- Compatible with single and king single bed sizes
- Rail can be installed on either side of the bed
- Easily adjustable in height to accommodate a range of mattress thicknesses, ensuring a secure fit
- Includes a convenient accessory pouch for storing personal items
- Provides sturdy and reliable support, making it easier to get in and out of bed

Code	Display Name	RRP	Box Quantity	Total Width	Total Height	Safe Working Load	Unit Weight
BEA013650	Affinity Grab Handle / Bed Rail	182	1	450mm	390 - 520mm	136kg	3.7kg