

## Folding Home Care Hoist 150kg - Sling separate



The Aspire 150F is a lightweight and portable folding patient lifter designed for home and community environments.



- Ergonomic Handset: The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Off Set Leg Design: Provides additional accessibility around chairs and furniture in the home
- Leg Function: Simple foot operation to open or close legs for lifting stability or accessibility
- Quick Release Yoke Attachment: One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety
- Build Quality: Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision
- Safety: Critical features include emergency stop button, emergency lowering and antientrapment technology, with both patient and caregiver safety at heart
- Intelligent Diagnostics: State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance

\*\*\* Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

**Code**

**Display Name**

H\_LSS390600 Folding Home Care Hoist 150kg - Sling separate

