

Folding Home Care Hoist 150kg - Sling separate



The Aspire 150F is a lightweight and portable folding patient lifter designed for home and community environments.



- **Ergonomic Handset:** The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- **Off Set Leg Design:** Provides additional accessibility around chairs and furniture in the home
- **Leg Function:** Simple foot operation to open or close legs for lifting stability or accessibility
- **Quick Release Yoke Attachment:** One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety
- **Build Quality:** Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision
- **Safety:** Critical features include emergency stop button, emergency lowering and antientrapment technology, with both patient and caregiver safety at heart
- **Intelligent Diagnostics:** State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance

*** Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

Code

Display Name

H_LSS390600 Folding Home Care Hoist 150kg - Sling separate

