

Liko Soft Original Highback Sling (Mod. 26)



- Soft pliable plate provides comfortable support when lifting from the recumbent position
- This sling offers a reclined seating position for patients with poor head and trunk stability
- Separated leg supports ensure that the patient does not sit on the sling after use
- Easily adaptable to the patient with comfortable full body support
- Can be used for floor transfers
- Available in Polyester

 Contract approved

ndis

Item Code:

05_123603164_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Warranty Information

- 1 Year



Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394560	Original Highback Sling - Reinforced Leg & Back Support - Polyester - MEDIUM SLIM	200 kg	Polyester (Loops)
LSS394570	Original Highback Sling - Reinforced Leg & Back Support - Polyester - MEDIUM	200 kg	Polyester (Loops)
LSS394580	Original Highback Sling - Reinforced Leg & Back Support - Polyester - LARGE	200 kg	Polyester (Loops)
LSS394590	Original Highback Sling - Reinforced Leg & Back Support - Polyester - XLARGE	500 kg	Polyester (Loops)
LSS394600	Original Highback Sling - Reinforced Leg & Back Support - Polyester - XXLARGE	500 kg	Polyester (Loops)