

## Liko Soft Original Highback Sling (Mod. 26)



- Soft pliable plate provides comfortable support when lifting from the recumbent position
- This sling offers a reclined seating position for patients with poor head and trunk stability
- Separated leg supports ensure that the patient does not sit on the sling after use
- Easily adaptable to the patient with comfortable full body support
- Can be used for floor transfers
- Available in Polyester



Contract approved

**ndis**

**Item Code:**

05\_123603164\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Warranty Information

- 1 Year



**Compliance**

AS/NZS ISO 10535:2011

| Code      | Display Name  | Safe Working Load | Material          |
|-----------|---|-------------------|-------------------|
| LSS394560 | Original Highback Sling - Reinforced Leg & Back Support - Polyester - MEDIUM SLIM | 200 kg            | Polyester (Loops) |
| LSS394570 | Original Highback Sling - Reinforced Leg & Back Support - Polyester - MEDIUM      | 200 kg            | Polyester (Loops) |
| LSS394580 | Original Highback Sling - Reinforced Leg & Back Support - Polyester - LARGE       | 200 kg            | Polyester (Loops) |
| LSS394590 | Original Highback Sling - Reinforced Leg & Back Support - Polyester - XLARGE      | 500 kg            | Polyester (Loops) |
| LSS394600 | Original Highback Sling - Reinforced Leg & Back Support - Polyester - XXLARGE     | 500 kg            | Polyester (Loops) |