

Milano Chair



- The Milano™ tilt in space chair is designed as a simple, adjustable and affordable care chair. It gives support, comfort and pressure management to patients who have low to medium support needs and those who still maintain a degree of mobility
- Tilt in space to allow for re-positioning for pressure care and comfort of the user. Elevating leg rest and recline are also standard
- Dartex multi stretch pressure management materials throughout all patient contact surfaces. This material is breathable with moisture permeability to ensure that the surface remains dry and comfortable for the user. This material can be wiped down for infection control and hygiene requirements
- Medium to high pressure management cushion, using memory foam. This cushion is removable for cleaning and if you wish to use a different type. An advantage of our memory foam cushions is they do not require maintenance or adjustment for each person, unlike a gel or air cushion. The cushion is also easily replaceable
- Adjustable seat width, foot plate height and seat depth allows the chair to fit the size and requirements of the user
- A light weight frame with directional castors to allow the chair to be easily pushed in a straight line or moved throughout a small room

Warranty Information

- 7 year frame
- 1 year upholstery

 Contract approved

Support @ Home

Function: Sitting furniture

Category: Under advice

Click [here](#) for more information

NDIS

Item Code:

05_180939182_0105_1_2

Click [here](#) for more information

Compliance

ISO13934-1 / ISO13937-1 / BS6807
/ BS7175

Code	Display Name	Total Length	Total Height	Safe Working Load	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Backrest Height
CHP197995	Pressure Relief Chair - Milano - Manual - Brown	-	-	-	-	-	-	-	-	-
CHP197997	Pressure Relief Chair - Milano - Electric - Taurus Cream	930 mm	1300 mm	159 kg	52 kg	Taurus Cream	400 - 500 mm	400 - 500 mm	540 mm	770 mm