


Liko Silhouette Sling (Mod. 22)



- Pliable sling that adapts to the body and requires minimal wheelchair space
- Recommended for lifting to sitting to/from close-fitting wheelchairs
- Provides a slightly reclined sitting posture suitable for patients with poor trunk stability
- Adjustable head support
- Suitable to be used during bathing
- Available in Mesh

Warranty Information

- 1 Year



ndis

Item Code:

05_123603164_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394450	Silhouette Mod 22 Sling - Mesh - XSMALL	200 kg	Mesh (Loops)
LSS394460	Silhouette Mod 22 Sling - Mesh - SMALL	-	-
LSS394470	Silhouette Mod 22 Sling - Mesh - MEDIUM	-	-
LSS394440	Silhouette Mod 22 Sling - Mesh - LARGE	-	-